

RTNETDEMO Example

Table of Contents

1 Symbol Reference 1

1.1 Alignment 1

1.1.1 Alignment::Alignment 1

1.1.2 Alignment::RunRTNet 1

1.2 AlignmentFactory 2

1.2.1 AlignmentFactory::AlignmentFactory 2

1.3 AcqContinuous 2

1.3.1 AcqContinuous::AcqContinuous 3

1.3.2 AcqContinuous::RunRTNet 3

1.4 AcqContinuousFactory 3

1.4.1 AcqContinuousFactory::AcqContinuousFactory 3

1.5 AcqContinuousBuffered 4

1.5.1 AcqContinuousBuffered::AcqContinuousBuffered 4

1.5.2 AcqContinuousBuffered::RunRTNet 4

1.6 AcqContinuousBufferedFactory 5

1.6.1 AcqContinuousBufferedFactory::AcqContinuousBufferedFactory 5

1.7 AcqSingleShot 5

1.7.1 AcqSingleShot::AcqSingleShot 6

1.7.2 AcqSingleShot::RunRTNet 6

1.8 AcqSingleShotFactory 6

1.8.1 AcqSingleShotFactory::AcqSingleShotFactory 7

1.9 CodaMode 7

1.9.1 CodaMode::CodaMode 7

1.9.2 CodaMode::RunRTNet 7

1.10 CodaModeFactory 8

1.10.1 CodaModeFactory::CodaModeFactory 8

1.11 StartSystem 8

1.11.1 StartSystem::RunRTNet 9

1.11.2 StartSystem::StartSystem 9

1.12 StartSystemFactory 9

1.12.1 StartSystemFactory::StartSystemFactory 9

1.13 RTNetDemoOptions 10

1.13.1 RTNetDemoOptions::CommandFile 10

1.13.2 RTNetDemoOptions::DataFile 10

1.13.3 RTNetDemoOptions::Parse 10

1.13.4 RTNetDemoOptions::RTNetDemoOptions 11

1.13.5 RTNetDemoOptions::Server 11

1.14 CommandRTNet 11

1.14.1 CommandRTNet::Client 12

| | | |
|-------------|--|-----------|
| 1.14.2 | CommandRTNet::Run | 12 |
| 1.14.3 | CommandRTNet::RunRTNet | 12 |
| 1.14.4 | friend class CommandFactoryRTNet | 12 |
| 1.15 | CommandFactoryRTNet | 12 |
| 1.15.1 | CommandFactoryRTNet::Build | 13 |
| 1.15.2 | CommandFactoryRTNet::CommandFactoryRTNet | 13 |
| 1.16 | CommandRTNetData | 13 |
| 1.16.1 | CommandRTNetData::DataOutput | 14 |
| 1.16.2 | friend class CommandFactoryRTNetData | 14 |
| 1.17 | CommandFactoryRTNetData | 14 |
| 1.17.1 | CommandFactoryRTNetData::Build | 15 |
| 1.17.2 | CommandFactoryRTNetData::CommandFactoryRTNetData | 15 |
| 1.18 | Framework | 15 |
| 1.18.1 | CommandLineOptions | 15 |
| 1.18.1.1 | CommandLineOptions::CommandLineOptions | 16 |
| 1.18.1.2 | CommandLineOptions::OptionRegister | 16 |
| 1.18.1.3 | CommandLineOptions::Parse | 16 |
| 1.18.1.4 | CommandLineOptions::RegisterOption | 16 |
| 1.18.2 | ResultLog | 17 |
| 1.18.2.1 | ResultLog::Log1 | 17 |
| 1.18.2.2 | ResultLog::Log2 | 17 |
| 1.18.2.3 | ResultLog::LogException | 18 |
| 1.18.2.4 | ResultLog::ResultLog | 18 |
| 1.18.3 | TracedException | 18 |
| 1.18.3.1 | TracedException::Report | 19 |
| 1.18.3.2 | TracedException::TracedException | 19 |
| 1.18.3.3 | TracedException::TracedException | 19 |
| 1.18.4 | STOP | 19 |
| 1.18.5 | STOP2 | 19 |
| 1.18.6 | Commands | 20 |
| 1.18.6.1 | Command | 20 |
| 1.18.6.1.1 | Command::~~Command | 21 |
| 1.18.6.1.2 | Command::Command | 21 |
| 1.18.6.1.3 | Command::Comment | 21 |
| 1.18.6.1.4 | Command::Comment | 21 |
| 1.18.6.1.5 | Command::Comment1 | 21 |
| 1.18.6.1.6 | Command::Comment1 | 22 |
| 1.18.6.1.7 | Command::GetParameter | 22 |
| 1.18.6.1.8 | Command::HaveRequiredParameters | 22 |
| 1.18.6.1.9 | Command::Name | 22 |
| 1.18.6.1.10 | Command::RegisterParameter | 22 |
| 1.18.6.1.11 | Command::Run | 22 |

| | |
|---|----|
| 1.18.6.1.12 friend class CommandFactory | 23 |
| 1.18.6.2 CommandFactory | 23 |
| 1.18.6.2.1 CommandFactory::Build | 23 |
| 1.18.6.2.2 CommandFactory::CommandFactory | 23 |
| 1.18.6.3 CommandList | 23 |
| 1.18.6.3.1 CommandList::~CommandList | 24 |
| 1.18.6.3.2 CommandList::BuildCommand | 24 |
| 1.18.6.3.3 CommandList::CommandList | 24 |
| 1.18.6.3.4 CommandList::RegisterCommand | 24 |
| 1.18.6.3.5 CommandList::RunAll | 25 |
| 1.18.6.4 Parameter | 25 |
| 1.18.6.4.1 Parameter::~Parameter | 25 |
| 1.18.6.4.2 Parameter::IsSet | 25 |
| 1.18.6.4.3 Parameter::Parameter | 25 |
| 1.18.6.4.4 Parameter::Parse | 26 |
| 1.18.6.5 ParameterInteger | 26 |
| 1.18.6.5.1 ParameterInteger::~ParameterInteger | 26 |
| 1.18.6.5.2 ParameterInteger::ParameterInteger | 26 |
| 1.18.6.5.3 ParameterInteger::Parse | 27 |
| 1.18.6.5.4 ParameterInteger::Value | 27 |
| 1.18.6.6 ParameterString | 27 |
| 1.18.6.6.1 ParameterString::~ParameterString | 27 |
| 1.18.6.6.2 ParameterString::ParameterString | 27 |
| 1.18.6.6.3 ParameterString::Parse | 28 |
| 1.18.6.6.4 ParameterString::Value | 28 |
| 1.18.6.7 COMMAND_STOP | 28 |
| 1.18.6.8 COMMAND_STOP2 | 28 |
| 1.18.7 IO | 29 |
| 1.18.7.1 DataAcquisition | 29 |
| 1.18.7.1.1 DataAcquisition::AddPacket | 29 |
| 1.18.7.1.2 DataAcquisition::DataAcquisition | 30 |
| 1.18.7.1.3 DataAcquisition::GetDeviceSampleRate | 30 |
| 1.18.7.1.4 DataAcquisition::Packets | 30 |
| 1.18.7.1.5 DataAcquisition::SetDeviceSampleRate | 30 |
| 1.18.7.2 DataPacket | 30 |
| 1.18.7.2.1 DataPacket::device | 31 |
| 1.18.7.2.2 DataPacket::page | 31 |
| 1.18.7.2.3 DataPacket::tick | 31 |
| 1.18.7.2.4 DataPacket::~ | 31 |
| 1.18.7.2.5 DataPacket::AddPoint | 31 |
| 1.18.7.2.6 DataPacket::DataPacket | 32 |
| 1.18.7.2.7 DataPacket::DataPacket | 32 |

| | |
|---|----|
| 1.18.7.2.8 DataPacket::Points | 32 |
| 1.18.7.3 DataPoint | 32 |
| 1.18.7.3.1 DataPoint::= | 33 |
| 1.18.7.3.2 DataPoint::AddIntensity | 33 |
| 1.18.7.3.3 DataPoint::AddResidual | 33 |
| 1.18.7.3.4 DataPoint::AddValue | 33 |
| 1.18.7.3.5 DataPoint::DataPoint | 33 |
| 1.18.7.3.6 DataPoint::DataPoint | 33 |
| 1.18.7.3.7 DataPoint::Intensity | 34 |
| 1.18.7.3.8 DataPoint::Occluded | 34 |
| 1.18.7.3.9 DataPoint::Residual | 34 |
| 1.18.7.3.10 DataPoint::SetOccluded | 34 |
| 1.18.7.3.11 DataPoint::Value | 34 |
| 1.18.7.4 DataReader | 34 |
| 1.18.7.4.1 DataReader::DataReader | 35 |
| 1.18.7.4.2 DataReader::LoadAcquisition | 35 |
| 1.18.7.4.3 DataReader::LocateObject | 35 |
| 1.18.7.4.4 DataReader::ParseFloat32 | 36 |
| 1.18.7.4.5 DataReader::ParseInt32 | 36 |
| 1.18.7.5 DataWriter | 36 |
| 1.18.7.5.1 DataWriter::AddAcquisitionPacket | 37 |
| 1.18.7.5.2 DataWriter::AddVersion | 37 |
| 1.18.7.5.3 DataWriter::BeginAcquisition | 37 |
| 1.18.7.5.4 DataWriter::BeginFile | 37 |
| 1.18.7.5.5 DataWriter::DataWriter | 37 |
| 1.18.7.5.6 DataWriter::EndAcquisition | 38 |
| 1.18.7.5.7 DataWriter::EndFile | 38 |
| 1.18.7.6 TextReader | 38 |
| 1.18.7.6.1 TextReader::NextNonSpace | 38 |
| 1.18.7.6.2 TextReader::ParseKeyValue | 39 |
| 1.18.7.6.3 TextReader::ParseQuotedString | 39 |
| 1.18.7.6.4 TextReader::ParseValue | 39 |
| 1.18.7.6.5 TextReader::SkipSpace | 40 |
| 1.18.7.6.6 TextReader::TextReader | 40 |
| 1.18.7.7 CommandListReader | 40 |
| 1.18.7.7.1 CommandListReader::CommandListReader | 40 |
| 1.18.7.7.2 CommandListReader::Read | 41 |
| 1.18.8 Tests | 41 |
| 1.18.8.1 SelfTest | 41 |
| 1.18.8.2 SelfTestCommandLine | 42 |
| 1.18.8.3 SelfTestDataReader | 42 |
| 1.18.8.4 SelfTestDataWriter | 42 |

| | |
|-----------------------------|----|
| 1.18.8.5 SelfTestScripts | 42 |
| 1.18.8.6 SelfTestModuleName | 43 |
| 1.18.8.7 SELF_TEST_ASSERT | 43 |
| 1.18.8.8 SELF_TEST_FAIL | 43 |

2 Index 44

RTNETDEMO Example

1 Symbol Reference

1.1 Alignment

Class Hierarchy

Command
 CommandRTNet
 Alignment

class Alignment : public CommandRTNet;

File



Alignment.h

Description

Demonstrate alignment

Members

Methods

| Method | Description |
|--|-----------------------------|
|  Alignment (see page 1) | Constructor used by factory |
|  RunRTNet (see page 1) | Implementation |

Legend

| | |
|---|---------|
|  | Method |
|  | virtual |

1.1.1 Alignment::Alignment

```
Alignment();
```

Description

Constructor used by factory

1.1.2 Alignment::RunRTNet

```
virtual void RunRTNet() throw(TracedException, codaRTNet::NetworkException, codaRTNet::DeviceStatusArray);
```

Description

Implementation

1.2 AlignmentFactory

Class Hierarchy

CommandFactory
 CommandFactoryRTNet
 AlignmentFactory

```
class AlignmentFactory : public CommandFactoryRTNet;
```

File


Alignment.h

Description

Factory for alignment commands

Members

Methods

| Method | Description |
|---|---|
|  AlignmentFactory (🔗 see page 2) | This is AlignmentFactory, a member of class AlignmentFactory. |

Legend

| | |
|---|--------|
|  | Method |
|---|--------|

1.2.1 AlignmentFactory::AlignmentFactory

```
AlignmentFactory(ResultLog& _results, codaRTNet::RTNetClient& _client);
```

Description

This is AlignmentFactory, a member of class AlignmentFactory.

1.3 AcqContinuous

Class Hierarchy

Command
 CommandRTNet
 CommandRTNetData
 AcqContinuous

```
class AcqContinuous : public CommandRTNetData;
```

File



AcqContinuous.h

Description

Demonstrate use of continuous acquisition

Members

Methods

| Method | Description |
|--|-----------------------------|
|  AcqContinuous (🔗 see page 3) | Constructor used by factory |
|  RunRTNet (🔗 see page 3) | Implementation |

Legend

| | |
|---|---------|
|  | Method |
|  | virtual |

1.3.1 AcqContinuous::AcqContinuous

```
AcqContinuous();
```

Description

Constructor used by factory

1.3.2 AcqContinuous::RunRTNet

```
virtual void RunRTNet() throw(TracedException, codaRTNet::NetworkException,
codaRTNet::DeviceStatusArray);
```

Description

Implementation

1.4 AcqContinuousFactory

Class Hierarchy

```
CommandFactory
  CommandFactoryRTNet
    CommandFactoryRTNetData
      AcqContinuousFactory
```

```
class AcqContinuousFactory : public CommandFactoryRTNetData;
```

File


AcqContinuous.h

Description

Factory for continuous acquisition commands

Members

Methods

| Method | Description |
|---|---|
|  AcqContinuousFactory (see page 3) | This is AcqContinuousFactory, a member of class AcqContinuousFactory. |

Legend

| | |
|---|--------|
|  | Method |
|---|--------|

1.4.1 AcqContinuousFactory::AcqContinuousFactory

```
AcqContinuousFactory(ResultLog& _results, codaRTNet::RTNetClient& _client, DataWriter&
_dataoutput);
```

Description

This is AcqContinuousFactory, a member of class AcqContinuousFactory.

1.5 AcqContinuousBuffered

Class Hierarchy

```
Command
  CommandRTNet
    CommandRTNetData
      AcqContinuousBuffered

class AcqContinuousBuffered : public CommandRTNetData;
```

File

AcqContinuousBuffered.h

Description

Demonstrate use of continuous buffered acquisition

Members

Methods

| Method | Description |
|---|-----------------------------|
|  AcqContinuousBuffered (see page 4) | Constructor used by factory |
|   RunRTNet (see page 4) | Implementation |

Legend

| | |
|---|---------|
|  | Method |
|  | virtual |

1.5.1 AcqContinuousBuffered::AcqContinuousBuffered

```
AcqContinuousBuffered();
```

Description

Constructor used by factory

1.5.2 AcqContinuousBuffered::RunRTNet

```
virtual void RunRTNet() throw(TracedException, codaRTNet::NetworkException,
codaRTNet::DeviceStatusArray);
```

Description

Implementation

1.6 AcqContinuousBufferedFactory

Class Hierarchy

```
CommandFactory
├── CommandFactoryRTNet
│   ├── CommandFactoryRTNetData
│   │   └── AcqContinuousBufferedFactory
│
└── class AcqContinuousBufferedFactory : public CommandFactoryRTNetData;
```

File


AcqContinuousBuffered.h

Description

Factory for continuous buffered acquisition commands

Members

Methods

| Method | Description |
|---|---|
|  AcqContinuousBufferedFactory (see page 5) | This is AcqContinuousBufferedFactory, a member of class AcqContinuousBufferedFactory. |

Legend

| | |
|---|--------|
|  | Method |
|---|--------|

1.6.1 AcqContinuousBufferedFactory::AcqContinuousBufferedFactory

```
AcqContinuousBufferedFactory(ResultLog& _results, codaRTNet::RTNetClient& _client,
    DataWriter& _dataoutput);
```

Description

This is AcqContinuousBufferedFactory, a member of class AcqContinuousBufferedFactory.

1.7 AcqSingleShot

Class Hierarchy

```
Command
├── CommandRTNet
│   ├── CommandRTNetData
│   │   └── AcqSingleShot
│
└── class AcqSingleShot : public CommandRTNetData;
```

File



AcqSingleShot.h

Description

Demonstrate single-shot acquisition command

Members

Methods

| Method | Description |
|--|-----------------------------|
|  AcqSingleShot (see page 6) | Constructor used by factory |
|  RunRTNet (see page 6) | Implementation |

Legend

| | |
|---|---------|
|  | Method |
|  | virtual |

1.7.1 AcqSingleShot::AcqSingleShot

```
AcqSingleShot();
```

Description

Constructor used by factory

1.7.2 AcqSingleShot::RunRTNet

```
virtual void RunRTNet() throw(TracedException, codaRTNet::NetworkException, codaRTNet::DeviceStatusArray);
```

Description

Implementation

1.8 AcqSingleShotFactory

Class Hierarchy

```
CommandFactory
  CommandFactoryRTNet
    CommandFactoryRTNetData
      AcqSingleShotFactory
```

```
class AcqSingleShotFactory : public CommandFactoryRTNetData;
```

File

AcqSingleShot.h

Description

Factory for single-shot acquisition commands

Members

Methods

| Method | Description |
|---|---|
|  AcqSingleShotFactory (see page 7) | This is AcqSingleShotFactory, a member of class AcqSingleShotFactory. |

Legend

| | |
|---|--------|
|  | Method |
|---|--------|

1.8.1 AcqSingleShotFactory::AcqSingleShotFactory

```
AcqSingleShotFactory(ResultLog& _results, codaRTNet::RTNetClient& _client, DataWriter&
_dataoutput);
```

Description
This is AcqSingleShotFactory, a member of class AcqSingleShotFactory.

1.9 CodaMode

Class Hierarchy

Command

CommandRTNet

CodaMode

```
class CodaMode : public CommandRTNet;
```

File



CodaMode.h

Description

Demonstrate setting of CODA mode

Members

Methods

| Method | Description |
|---|-----------------------------|
|  CodaMode (see page 7) | Constructor used by factory |
|  RunRTNet (see page 7) | Implementation |

Legend

| | |
|---|---------|
|  | Method |
|  | virtual |

1.9.1 CodaMode::CodaMode

```
CodaMode();
```

Description
Constructor used by factory

1.9.2 CodaMode::RunRTNet

```
virtual void RunRTNet() throw(TracedException, codaRTNet::NetworkException,
codaRTNet::DeviceStatusArray);
```

Description
Implementation

1.10 CodaModeFactory

Class Hierarchy

CommandFactory
CommandFactoryRTNet
CodaModeFactory

```
class CodaModeFactory : public CommandFactoryRTNet;
```

File


CodaMode.h

Description

Factory for code mode commands

Members

Methods

| Method | Description |
|--|---|
|  CodaModeFactory (🔗 see page 8) | This is CodaModeFactory, a member of class CodaModeFactory. |

Legend

| | |
|---|--------|
|  | Method |
|---|--------|

1.10.1 CodaModeFactory::CodaModeFactory

```
CodaModeFactory(ResultLog& _results, codaRTNet::RTNetClient& _client);
```

Description

This is CodaModeFactory, a member of class CodaModeFactory.

1.11 StartSystem

Class Hierarchy

Command
CommandRTNet
StartSystem

```
class StartSystem : public CommandRTNet;
```

File




StartSystem.h

Description

Demonstrate system start-up

Members

Methods

| Method | Description |
|---|-----------------------------|
|   RunRTNet (🔗 see page 9) | Implementation |
|  StartSystem (🔗 see page 9) | Constructor used by factory |

Legend

| | |
|---|---------|
|  | Method |
|  | virtual |

1.11.1 StartSystem::RunRTNet

```
virtual void RunRTNet() throw(TracedException, codaRTNet::NetworkException,
codaRTNet::DeviceStatusArray);
```

Description

Implementation

1.11.2 StartSystem::StartSystem

```
StartSystem();
```

Description

Constructor used by factory

1.12 StartSystemFactory

Class Hierarchy

```
CommandFactory
  CommandFactoryRTNet
    StartSystemFactory
```

```
class StartSystemFactory : public CommandFactoryRTNet;
```

File


StartSystem.h

Description

Factory to make start-system commands

Members

Methods

| Method | Description |
|---|---|
|  StartSystemFactory (see page 9) | This is StartSystemFactory, a member of class StartSystemFactory. |

Legend

| | |
|---|--------|
|  | Method |
|---|--------|

1.12.1 StartSystemFactory::StartSystemFactory

```
StartSystemFactory(ResultLog& _results, codaRTNet::RTNetClient& _client);
```

Description

This is StartSystemFactory, a member of class StartSystemFactory.

1.13 RTNetDemoOptions

Class Hierarchy

CommandLineOptions
RTNetDemoOptions

```
class RTNetDemoOptions : protected CommandLineOptions;
```

File

RTNetDemoOptions.h

Description

Command (see page 20) line options class specific to our RTNetDemo program

Members

Methods

| Method | Description |
|----------------------------------|---|
| ◆ CommandFile (see page 10) | Structured text file containing commands |
| ◆ DataFile (see page 10) | File to write data output |
| ◆ Parse (see page 10) | Call CommandLineOptions::Parse (see page 16) and log the resulting set of options |
| ◆ RTNetDemoOptions (see page 11) | This is RTNetDemoOptions, a member of class RTNetDemoOptions. |
| ◆ Server (see page 11) | RT Net server ip address as string |

Legend

| | |
|---|--------|
| ◆ | Method |
|---|--------|

1.13.1 RTNetDemoOptions::CommandFile

```
const std::string& CommandFile() const;
```

Description

Structured text file containing commands

1.13.2 RTNetDemoOptions::DataFile

```
const std::string& DataFile() const;
```

Description

File to write data output

1.13.3 RTNetDemoOptions::Parse

```
void Parse(ResultLog& results, int argc, char* argv[]) throw(TracedException);
```

Parameters

| Parameters | Description |
|--------------------|-------------------------|
| ResultLog& results | Log to write options to |
| int argc | The argc from main |

| | |
|--------------|--------------------|
| char* argv[] | The argv from main |
|--------------|--------------------|

Description

Call CommandLineOptions::Parse (see page 16) and log the resulting set of options

Exceptions

TracedException (see page 18) if there was a problem parsing options in base class

1.13.4 RTNetDemoOptions::RTNetDemoOptions

```
RTNetDemoOptions();
```

Description

This is RTNetDemoOptions, a member of class RTNetDemoOptions.

1.13.5 RTNetDemoOptions::Server

```
const std::string& Server() const;
```

Description

RT Net server ip address as string

1.14 CommandRTNet

Class Hierarchy

```
Command
  CommandRTNet
```

```
class CommandRTNet : public Command;
```

File

CommandRTNet.h

Description

Extension of abstract Command (see page 20) class to allow global RT Net SDK client to be set from factory and RT Net SDK exceptions to be caught

Members**Methods**

| Method | Description |
|------------------------|--|
| Client (see page 12) | Get client which was set by factory just after construction |
| Run (see page 12) | Implementation of abstract Run method to use new abstract RunRTNet (see page 12) method and catch RTNet exceptions |
| RunRTNet (see page 12) | Abstract method to be used by derived classes |

Friends

| Friend | Description |
|---|--|
| class CommandFactoryRTNet (see page 12) | This is friend friend class CommandFactoryRTNet. |

Legend

| | |
|---|----------|
|  | Method |
|  | virtual |
|  | abstract |

1.14.1 CommandRTNet::Client

```
codaRTNet::RTNetClient& Client();
```

Description

Get client which was set by factory just after construction

1.14.2 CommandRTNet::Run

```
virtual void Run() throw(TracedException);
```

Description

Implementation of abstract Run method to use new abstract RunRTNet (see page 12) method and catch RTNet exceptions

1.14.3 CommandRTNet::RunRTNet

```
virtual void RunRTNet() = 0 throw(codaRTNet::NetworkException,
codaRTNet::DeviceStatusArray);
```

Description

Abstract method to be used by derived classes

1.14.4 friend class CommandFactoryRTNet

```
friend class CommandFactoryRTNet;
```

Description

This is friend friend class CommandFactoryRTNet.

1.15 CommandFactoryRTNet

Class Hierarchy

```
CommandFactory
  CommandFactoryRTNet
```

```
class CommandFactoryRTNet : public CommandFactory;
```



File

CommandRTNet.h

Description

Extension of command factory to allow RT Net SDK client to be set in the commands we build

Members**Methods**

| Method | Description |
|---|--|
|  Build (see page 13) | Build new command and set RT Net SDK client |
|  CommandFactoryRTNet (see page 13) | Construct knowing the log and client which we will give to the commands we build |

Legend

| | |
|---|---------|
|  | Method |
|  | virtual |

1.15.1 CommandFactoryRTNet::Build

```
virtual CommandRTNet* Build(const std::string& name);
```

Parameters

| Parameters | Description |
|-------------------------|--------------------------------|
| const std::string& name | Name to use for command module |

Description

Build new command and set RT Net SDK client

1.15.2 CommandFactoryRTNet::CommandFactoryRTNet

```
CommandFactoryRTNet(ResultLog& _results, codaRTNet::RTNetClient& _client);
```

Parameters

| Parameters | Description |
|---------------------------------|---------------------------------------|
| ResultLog& _results | Results log to give to commands |
| codaRTNet::RTNetClient& _client | RT Net SDK client to give to commands |

Description

Construct knowing the log and client which we will give to the commands we build

1.16 CommandRTNetData

Class Hierarchy

```

Command
├── CommandRTNet
│   └── CommandRTNetData
```

```
class CommandRTNetData : public CommandRTNet;
```


File

CommandRTNetData.h

Description

Extend abstract RTNet command class also to allow a data output writer to be set by factory

Members**Methods**

| Method | Description |
|--|---|
|  DataOutput (see page 14) | Get data writer as set by factory just after construction |

Friends

| Friend | Description |
|---|--|
| class CommandFactoryRTNetData (see page 14) | This is friend friend class CommandFactoryRTNetData. |

Legend

| | |
|---|--------|
|  | Method |
|---|--------|

1.16.1 CommandRTNetData::DataOutput

```
DataWriter& DataOutput();
```

Description

Get data writer as set by factory just after construction

1.16.2 friend class CommandFactoryRTNetData

```
friend class CommandFactoryRTNetData;
```

Description

This is friend friend class CommandFactoryRTNetData.

1.17 CommandFactoryRTNetData

Class Hierarchy

```
CommandFactory
  CommandFactoryRTNet
    CommandFactoryRTNetData
```

```
class CommandFactoryRTNetData : public CommandFactoryRTNet;
```




File

CommandRTNetData.h

Description

Extend general RT Net command factory to allow data output writer to be set in the commands we build

Members**Methods**

| Method | Description |
|---|--|
|   Build (see page 15) | Build new command and set data writer |
|  CommandFactoryRTNetData (see page 15) | Construct knowing the log, client, and data output which we will give to the commands we build |

Legend

| | |
|---|---------|
|  | Method |
|  | virtual |

1.17.1 CommandFactoryRTNetData::Build

```
virtual CommandRTNetData* Build(const std::string& name);
```

Parameters

| Parameters | Description |
|-------------------------|--------------------------------|
| const std::string& name | Name to use for command module |

Description

Build new command and set data writer

1.17.2

CommandFactoryRTNetData::CommandFactoryRTNetData

```
CommandFactoryRTNetData(ResultLog& _results, codaRTNet::RTNetClient& _client, DataWriter& _dataoutput);
```

Parameters

| Parameters | Description |
|---------------------------------|---------------------------------------|
| ResultLog& _results | Results log to give to commands |
| codaRTNet::RTNetClient& _client | RT Net SDK client to give to commands |

Description

Construct knowing the log, client, and data output which we will give to the commands we build

1.18 Framework

Framework classes used by the RTNetDemo example

Classes

| Class | Description |
|----------------------------------|--|
| CommandLineOptions (see page 15) | Register command line options to look for and parse the command line accordingly. |
| ResultLog (see page 17) | Produce a formatted log |
| TracedException (see page 18) | An exception class which can log the source file and source code line at which the exception was thrown, plus module name and error message. |

Macros

| Macro | Description |
|---------------------|---|
| STOP (see page 19) | Throw a TracedException (see page 18) tagged with this line of source code and source file name. |
| STOP2 (see page 19) | Throw a TracedException (see page 18) tagged with this line of source code and source file name, and two message strings which are concatenated and separated by a colon. |

1.18.1 CommandLineOptions

Class Hierarchy

CommandLineOptions

```
class CommandLineOptions;
```





File

CommandLineOptions.h

Description

Register command line options to look for and parse the command line accordingly.

Members**Methods**

| Method | Description |
|--|---|
|  CommandLineOptions (see page 16) | Constructor (constructs with no options) |
|  OptionRegister (see page 16) | Get the array of all options registered |
|  Parse (see page 16) | Parse a command line |
|  RegisterOption (see page 16) | Register a possible option with the default value to use if not found |

Legend

| | |
|---|--------|
|  | Method |
|---|--------|

1.18.1.1 CommandLineOptions::CommandLineOptions

```
CommandLineOptions();
```

Description

Constructor (constructs with no options)

1.18.1.2 CommandLineOptions::OptionRegister

```
const std::map<std::string, std::string*>& OptionRegister() const;
```

Returns

Array of options

Description

Get the array of all options registered

1.18.1.3 CommandLineOptions::Parse

```
void Parse(int argc, char* argv[]) throw(TracedException);
```

Parameters

| Parameters | Description |
|--------------|--------------------|
| int argc | The argc from main |
| char* argv[] | The argv from main |

Description

Parse a command line

Exceptions

TracedException ([see page 18](#))

1.18.1.4 CommandLineOptions::RegisterOption

```
void RegisterOption(const char* name, std::string& option, const char* defaultvalue);
```

Parameters

| Parameters | Description |
|---------------------------------------|--|
| <code>const char* name</code> | Name to find on command line (the name on the command line will be prefixed with --) |
| <code>std::string& option</code> | A string object to fill with the value found for this option |
| <code>const char* defaultvalue</code> | The default value to put in option if no value found |

Description

Register a possible option with the default value to use if not found

1.18.2 ResultLog

Class Hierarchy

ResultLog

```
class ResultLog;
```





File

ResultLog.h

Description

Produce a formatted log

Members**Methods**

| Method | Description |
|--|--|
|  Log1 (see page 17) | Log a message |
|  Log2 (see page 17) | Concatenate two messages and log them |
|  LogException (see page 18) | Log the given exception which includes source file name and line of code |
|  ResultLog (see page 18) | Construct log to specified stream |

Legend

| | |
|---|--------|
|  | Method |
|---|--------|

1.18.2.1 ResultLog::Log1

```
void Log1(const std::string& module, const std::string& message);
```

Parameters

| Parameters | Description |
|---|---|
| <code>const std::string& module</code> | Module name to use in prefix to message |
| <code>const std::string& message</code> | Message to write |

Description

Log a message

1.18.2.2 ResultLog::Log2

```
void Log2(const std::string& module, const std::string& message1, const std::string& message2);
```

Parameters

| Parameters | Description |
|--|---|
| <code>const std::string& module</code> | Module name to use in prefix to message |
| <code>const std::string& message1</code> | First message to write, will be followed by colon and space |

| | |
|--|-------------------------|
| <code>const std::string& message2</code> | Second message to write |
|--|-------------------------|

Description

Concatenate two messages and log them

1.18.2.3 ResultLog::LogException

```
void LogException(const std::string& module, const TracedException& e);
```

Parameters

| Parameters | Description |
|--|--|
| <code>const std::string& module</code> | Module to use in prefix (should be the logging module, not necessarily the module which threw the exception) |
| <code>const TracedException& e</code> | Exception to log |

Description

Log the given exception which includes source file name and line of code

1.18.2.4 ResultLog::ResultLog

```
ResultLog(std::ostream& _diag);
```

Parameters

| Parameters | Description |
|--------------------------------------|------------------|
| <code>std::ostream& _diag</code> | Stream to log to |

Description

Construct log to specified stream

1.18.3 TracedException

An exception class which can log the source file and source code line at which the exception was thrown, plus module name and error message.

Class Hierarchy

TracedException

```
class TracedException;
```




File

TracedException.h

Description

This class is not to be thrown directly but rather via the STOP (see page 19) and STOP2 (see page 19) macros which use the C preprocessor to log the source file and source line at which they are used.

Members**Methods**

| Method | Description |
|---|--|
|  Report (see page 19) | Print this exception to the given stream |
|  TracedException (see page 19) | Copy constructor |
|  TracedException (see page 19) | Constructor to be called from STOP (see page 19) and STOP2 (see page 19) |

Legend

| | |
|---|--------|
|  | Method |
|---|--------|

1.18.3.1 TracedException::Report

```
void Report(std::ostream& report) const;
```

Parameters

| Parameters | Description |
|----------------------|------------------------|
| std::ostream& report | The stream to print to |

Description

Print this exception to the given stream

1.18.3.2 TracedException::TracedException

```
TracedException(const TracedException& src);
```

Parameters

| Parameters | Description |
|----------------------------|------------------|
| const TracedException& src | Instance to copy |

Description

Copy constructor

1.18.3.3 TracedException::TracedException

```
TracedException(const std::string& _module, const char* _sourcefile, unsigned long
_sourcecline, const std::string& _message);
```

Description

Constructor to be called from STOP (see page 19) and STOP2 (see page 19)

1.18.4 STOP

```
#define STOP(_module, _message) throw TracedException(_module, __FILE__, __LINE__, _message)
```

File

TracedException.h

Parameters

| Parameters | Description |
|------------|---------------------------------------|
| _module | Module name to put in the exception |
| _message | Error message to put in the exception |

Description

Throw a TracedException (see page 18) tagged with this line of source code and source file name.

1.18.5 STOP2

```
#define STOP2(_module, _message1, _value) STOP(_module, std::string(_message1) +
std::string(": ") + _value)
```

File

TracedException.h

Parameters

| Parameters | Description |
|------------|--|
| _module | Module name to put in the exception |
| _message1 | Error message to put in the exception |
| _value | Value to append to error message following a colon |

Description

Throw a TracedException (see page 18) tagged with this line of source code and source file name, and two message strings which are concatenated and separated by a colon.

1.18.6 Commands

Classes to allow a list of commands to be created and run.

Classes

| Class | Description |
|--------------------------------|--|
| Command (see page 20) | The purpose of this class is to have an abstract way to describe commands such that a list of them can be read from a file and executed. Derive from this class to add scriptable commands to the program. Parameters can be added to the command by having member variables derived from Parameter (see page 25) which are registered using the RegisterParameter (see page 22) method. |
| CommandFactory (see page 23) | A factory object for building commands of a specified type Derive from this to implement factories for specific commands |
| CommandList (see page 23) | A list of command factories and the ability to manufacture commands by name and keep them in a local list |
| Parameter (see page 25) | Base class for parameters which we can register with our command objects |
| ParameterInteger (see page 26) | Integer command parameter |
| ParameterString (see page 27) | Text string command parameter |

Macros

| Macro | Description |
|-----------------------------|---|
| COMMAND_STOP (see page 28) | Call the STOP (see page 19) macro to throw a TracedException (see page 18) to log this source file and line plus this command's name |
| COMMAND_STOP2 (see page 28) | Call the STOP2 (see page 19) macro to throw a TracedException (see page 18) to log this source file and line plus this command's name. Also like STOP2 (see page 19) this will concatenate message and value. |

1.18.6.1 Command

The purpose of this class is to have an abstract way to describe commands such that a list of them can be read from a file and executed. Derive from this class to add scriptable commands to the program. Parameters can be added to the command by having member variables derived from Parameter (see page 25) which are registered using the RegisterParameter (see page 22) method.

Class Hierarchy

Command












class Command;**File**

Command.h

Description

An abstract command which can be executed by our program.



Members**Methods**

| Method | Description |
|---|--|
|  ~Command (see page 21) | Virtual destructor |
|  Command (see page 21) | Protected empty constructor to be used by the CommandFactory (see page 23) |
|  Comment (see page 21) | Add message to log using this command's name as module name |
|  Comment (see page 21) | Add contents of string stream to log using this command's name as module name |
|  Comment1 (see page 21) | Add message and string to log using this command's name as module name |
|  Comment1 (see page 22) | Add message and number to log using this command's name as module name |
|  GetParameter (see page 22) | Get the wrapper object for a parameter |
|  HaveRequiredParameters (see page 22) | See if all required parameters are set |
|  Name (see page 22) | Name of this command |
|  RegisterParameter (see page 22) | Register a parameter name used by this command - useful for looking up parameters in script file |
|  Run (see page 22) | Abstract command method to be implemented by derived classes |

Friends

| Friend | Description |
|--|---|
| class CommandFactory (see page 23) | This is friend friend class CommandFactory. |

Legend

| | |
|---|-----------|
|  | Method |
|  | virtual |
|  | protected |
|  | abstract |

1.18.6.1.1 Command::~~Command

```
virtual ~Command();
```

Description

Virtual destructor

1.18.6.1.2 Command::Command

```
Command();
```

Description

Protected empty constructor to be used by the CommandFactory ([see page 23](#))

1.18.6.1.3 Command::Comment

```
void Comment(const std::string& message);
```

Description

Add message to log using this command's name as module name

1.18.6.1.4 Command::Comment

```
void Comment(std::ostream& formattedmessage);
```

Description

Add contents of string stream to log using this command's name as module name

1.18.6.1.5 Command::Comment1

```
void Comment1(const std::string& message, const std::string& value);
```

Description

Add message and string to log using this command's name as module name

1.18.6.1.6 Command::Comment1

```
void Comment1(const std::string& message, long value);
```

Description

Add message and number to log using this command's name as module name

1.18.6.1.7 Command::GetParameter

```
Parameter* GetParameter(const std::string& str);
```

Parameters

| Parameters | Description |
|------------------------|---------------------------|
| const std::string& str | name of parameter to find |

Returns

The parameter wrapper if found, NULL otherwise

Description

Get the wrapper object for a parameter

1.18.6.1.8 Command::HaveRequiredParameters

```
bool HaveRequiredParameters() const;
```

Returns

true if all parameters filled, false otherwise

Description

See if all required parameters are set

1.18.6.1.9 Command::Name

```
const std::string& Name() const;
```

Description

Name of this command

1.18.6.1.10 Command::RegisterParameter

```
void RegisterParameter(const std::string& parametername, Parameter& parameter);
```

Description

Register a parameter name used by this command - useful for looking up parameters in script file

1.18.6.1.11 Command::Run

```
virtual void Run() = 0 throw(TracedException);
```

Description

Abstract command method to be implemented by derived classes

1.18.6.1.12 friend class CommandFactory

```
friend class CommandFactory;
```

Description

This is friend friend class CommandFactory.

1.18.6.2 CommandFactory

Class Hierarchy

```
CommandFactory
class CommandFactory;
```

File



CommandFactory.h

Description


A factory object for building commands of a specified type Derive from this to implement factories for specific commands

Members

Methods

| Method | Description |
|---|----------------------------|
|  Build (see page 23) | Build a new command object |
|  CommandFactory (see page 23) | Construct |

Legend

| | |
|---|---------|
|  | Method |
|  | virtual |

1.18.6.2.1 CommandFactory::Build

```
virtual Command* Build(const std::string& name);
```

Description

Build a new command object

1.18.6.2.2 CommandFactory::CommandFactory

```
CommandFactory(ResultLog& _results);
```

Parameters

| Parameters | Description |
|---------------------|--|
| ResultLog& _results | Results log to give to the command objects we make |

Description

Construct

1.18.6.3 CommandList

Class Hierarchy

```
CommandList
class CommandList;
```







File

CommandList.h

Description

A list of command factories and the ability to manufacture commands by name and keep them in a local list

Members**Methods**

| Method | Description |
|--|---|
|   ~CommandList (see page 24) | Destructor |
|  BuildCommand (see page 24) | Look up factory with specified command name, make new command and store in list |
|  CommandList (see page 24) | Construct empty |
|  RegisterCommand (see page 24) | Register a command factory by name |
|  RunAll (see page 25) | Run all commands in local list (will have been made using BuildCommand (see page 24)) |

Legend

| | |
|---|---------|
|  | Method |
|  | virtual |

1.18.6.3.1 CommandList::~~CommandList

```
virtual ~CommandList();
```

Description

Destructor

1.18.6.3.2 CommandList::BuildCommand

```
Command* BuildCommand(const std::string& commandname);
```

Parameters

| Parameters | Description |
|--------------------------------|--------------|
| const std::string& commandname | Name to find |

Description

Look up factory with specified command name, make new command and store in list

1.18.6.3.3 CommandList::CommandList

```
CommandList();
```

Description

Construct empty

1.18.6.3.4 CommandList::RegisterCommand

```
void RegisterCommand(const std::string& name, CommandFactory* factory);
```

Parameters

| Parameters | Description |
|-------------------------|-----------------------|
| const std::string& name | Name to register |
| CommandFactory* factory | Corresponding factory |

Description

Register a command factory by name

1.18.6.3.5 CommandList::RunAll

```
void RunAll() const throw(TracedException);
```

Description

Run all commands in local list (will have been made using BuildCommand (see page 24))

Exceptions

TracedException (see page 18) if a command fails

1.18.6.4 Parameter

Class Hierarchy

Parameter

```
class Parameter;
```

File

Parameter.h

Description

Base class for parameters which we can register with our command objects

Members

Methods

| Method | Description |
|--|---|
|   ~Parameter (see page 25) | Destructor |
|  IsSet (see page 25) | See if successfully parsed |
|  Parameter (see page 25) | Construct empty / default value |
|   Parse (see page 26) | Parse from string read from structured text |

Legend

| | |
|---|----------|
|  | Method |
|  | virtual |
|  | abstract |

1.18.6.4.1 Parameter::~~Parameter

```
virtual ~Parameter();
```

Description

Destructor

1.18.6.4.2 Parameter::IsSet

```
bool IsSet() const;
```

Returns

true if successfully parsed

Description

See if successfully parsed

1.18.6.4.3 Parameter::Parameter

```
Parameter();
```

Description

Construct empty / default value

1.18.6.4 Parameter::Parse

```
virtual void Parse(const std::string& str) = 0;
```

Parameters

| Parameters | Description |
|------------------------|---------------|
| const std::string& str | Text to parse |

Description

Parse from string read from structured text

1.18.6.5 ParameterInteger

Class Hierarchy

```
Parameter
  ParameterInteger

class ParameterInteger : public Parameter;
```

File





ParameterInteger.h

Description

Integer command parameter

Members

Methods

| Method | Description |
|---|----------------------------------|
|  ~ParameterInteger (see page 26) | Destructor |
|  ParameterInteger (see page 26) | Construct and initialise to zero |
|  Parse (see page 27) | Parse string to integer |
|  Value (see page 27) | Get value of integer |

Legend

| | |
|---|---------|
|  | Method |
|  | virtual |

1.18.6.5.1 ParameterInteger::~~ParameterInteger

```
virtual ~ParameterInteger();
```

Description

Destructor

1.18.6.5.2 ParameterInteger::ParameterInteger

```
ParameterInteger();
```

Description

Construct and initialise to zero

1.18.6.5.3 ParameterInteger::Parse

```
virtual void Parse(const std::string& str);
```

Description
Parse string to integer

1.18.6.5.4 ParameterInteger::Value

```
const long Value() const;
```

Returns
value

Description
Get value of integer

1.18.6.6 ParameterString

Class Hierarchy
Parameter
 ParameterString






```
class ParameterString : public Parameter;
```

File
ParameterString.h


Description
Text string command parameter

Members

Methods

| Method | Description |
|--|---------------------------------|
|  ~ParameterString (see page 27) | Destructor |
|  ParameterString (see page 27) | Construct empty |
|  Parse (see page 28) | Parse from string (just copies) |
|  Value (see page 28) | Get value |

Legend

| | |
|---|---------|
|  | Method |
|  | virtual |

1.18.6.6.1 ParameterString::~~ParameterString

```
virtual ~ParameterString();
```

Description
Destructor

1.18.6.6.2 ParameterString::ParameterString

```
ParameterString();
```

Description
Construct empty

1.18.6.6.3 ParameterString::Parse

```
virtual void Parse(const std::string& str);
```

Parameters

| Parameters | Description |
|------------------------|----------------|
| const std::string& str | String to copy |

Description

Parse from string (just copies)

1.18.6.6.4 ParameterString::Value

```
const std::string& Value() const;
```

Returns

text value

Description

Get value

1.18.6.7 COMMAND_STOP

```
#define COMMAND_STOP(_message) STOP(Name(), _message)
```

File

Command.h

Parameters

| Parameters | Description |
|------------|----------------------|
| _message | Error message to log |

Description

Call the STOP (see page 19) macro to throw a TracedException (see page 18) to log this source file and line plus this command's name

Exceptions

TracedException (see page 18)

1.18.6.8 COMMAND_STOP2

```
#define COMMAND_STOP2(_message, _value) STOP2(Name(), _message, _value)
```

File

Command.h

Parameters

| Parameters | Description |
|------------|----------------------|
| _message | Error message to log |

Description

Call the STOP2 (see page 19) macro to throw a TracedException (see page 18) to log this source file and line plus this command's name. Also like STOP2 (see page 19) this will concatenate message and value.

Exceptions

TracedException (see page 18)

1.18.7 IO

Input and output to text files

Classes

| Class | Description |
|---------------------------------|---|
| DataAcquisition (see page 29) | Data acquisition as read from structured text |
| DataPacket (see page 30) | Data packet as read from structured text corresponding to RTNetworkPacket |
| DataPoint (see page 32) | Data point (element of packet) as read from structured text |
| DataReader (see page 34) | Read data from structured text as output by DataWriter (see page 36) |
| DataWriter (see page 36) | Write data to a structured text file |
| TextReader (see page 38) | Read structured text files |
| CommandListReader (see page 40) | Read (see page 41) list of commands from structured text file |

1.18.7.1 DataAcquisition

Class Hierarchy

```
DataAcquisition
class DataAcquisition;
```

File

DataReader.h

Description

Data acquisition as read from structured text

Members

Methods

| Method | Description |
|-------------------------------------|--|
| ◆ AddPacket (see page 29) | Append a packet corresponding to RTNetworkPacket information |
| ◆ DataAcquisition (see page 30) | Construct empty |
| ◆ GetDeviceSampleRate (see page 30) | Retrieve sample rate info |
| ◆ Packets (see page 30) | Get array of packets |
| ◆ SetDeviceSampleRate (see page 30) | Append device sample-rate settings |

Legend

| | |
|---|--------|
| ◆ | Method |
|---|--------|

1.18.7.1.1 DataAcquisition::AddPacket

```
void AddPacket(DataPacket& packet);
```

Parameters

| Parameters | Description |
|--------------------|------------------|
| DataPacket& packet | Packet to append |

Description

Append a packet corresponding to RTNetworkPacket information

1.18.7.1.2 DataAcquisition::DataAcquisition

```
DataAcquisition();
```

Description

Construct empty

1.18.7.1.3 DataAcquisition::GetDeviceSampleRate

```
bool GetDeviceSampleRate(long device, float& rate);
```

Parameters

| Parameters | Description |
|-------------|------------------------------|
| long device | Device whose rate to look up |
| float& rate | Filled with rate if found |

Returns

true if found, false otherwise

Description

Retrieve sample rate info

1.18.7.1.4 DataAcquisition::Packets

```
const std::vector<DataPacket>& Packets() const;
```

Description

Get array of packets

1.18.7.1.5 DataAcquisition::SetDeviceSampleRate

```
void SetDeviceSampleRate(long device, float rate);
```

Parameters

| Parameters | Description |
|-------------|---------------|
| long device | Device to set |
| float rate | Rate in Hz |

Description

Append device sample-rate settings

1.18.7.2 DataPacket

Class Hierarchy

DataPacket

```
class DataPacket;
```

File




DataReader.h

Description






Data packet as read from structured text corresponding to RTNetworkPacket

Members

Data Members

| Data Member | Description |
|--|----------------------------------|
|  device (see page 31) | Device ID from RTNetworkPacket |
|  page (see page 31) | Device page from RTNetworkPacket |
|  tick (see page 31) | Device tick from RTNetworkPacket |

Methods

| Method | Description |
|--|--|
|  = (see page 31) | Copy operator |
|  AddPoint (see page 31) | Append a point |
|  DataPacket (see page 32) | Construct empty |
|  DataPacket (see page 32) | Copy constructor - calls copy operator |
|  Points (see page 32) | Get points array |

Legend

| | |
|---|-------------|
|  | Data Member |
|  | Method |

1.18.7.2.1 DataPacket::device

```
long device;
```

Description

Device ID from RTNetworkPacket

1.18.7.2.2 DataPacket::page

```
long page;
```

Description

Device page from RTNetworkPacket

1.18.7.2.3 DataPacket::tick

```
long tick;
```

Description

Device tick from RTNetworkPacket

1.18.7.2.4 DataPacket::=

```
const DataPacket& operator =(const DataPacket& src);
```

Description

Copy operator

1.18.7.2.5 DataPacket::AddPoint

```
void AddPoint(const DataPoint& p);
```

Parameters

| Parameters | Description |
|--------------------|-----------------|
| const DataPoint& p | Point to append |

Description

Append a point

1.18.7.2.6 DataPacket::DataPacket

DataPacket () ;

Description

Construct empty

1.18.7.2.7 DataPacket::DataPacket

DataPacket (const DataPacket& src) ;

Description

Copy constructor - calls copy operator

1.18.7.2.8 DataPacket::Points

const std::vector<DataPoint>& Points() const ;

Description

Get points array

1.18.7.3 DataPoint

Class Hierarchy

DataPoint

class DataPoint ;

File












DataReader.h

Description

Data point (element of packet) as read from structured text

Members

Methods

| Method | Description |
|--|---|
|  = (see page 33) | Copy operator |
|  AddIntensity (see page 33) | Append intensity value |
|  AddResidual (see page 33) | Append residual value |
|  AddValue (see page 33) | Append coordinate value |
|  DataPoint (see page 33) | Construct empty |
|  DataPoint (see page 33) | Copy constructor - calls copy operator=() |
|  Intensity (see page 34) | Get intensity array |
|  Occluded (see page 34) | Get occlusion flag |
|  Residual (see page 34) | Get residual array |
|  SetOccluded (see page 34) | Set occlusion flag |
|  Value (see page 34) | Get value array |

Legend

| | |
|---|--------|
|  | Method |
|---|--------|

1.18.7.3.1 DataPoint::=

```
const DataPoint& operator =(const DataPoint& src);
```

Description

Copy operator

1.18.7.3.2 DataPoint::AddIntensity

```
void AddIntensity(long i);
```

Parameters

| Parameters | Description |
|------------|-------------|
| long i | intensity |

Description

Append intensity value

1.18.7.3.3 DataPoint::AddResidual

```
void AddResidual(long r);
```

Parameters

| Parameters | Description |
|------------|-------------|
| long r | residual |

Description

Append residual value

1.18.7.3.4 DataPoint::AddValue

```
void AddValue(float x);
```

Parameters

| Parameters | Description |
|------------|---------------------|
| float x | Value (see page 34) |

Description

Append coordinate value

1.18.7.3.5 DataPoint::DataPoint

```
DataPoint();
```

Description

Construct empty

1.18.7.3.6 DataPoint::DataPoint

```
DataPoint(const DataPoint& src);
```

Description

Copy constructor - calls copy operator=()

1.18.7.3.7 DataPoint::Intensity

```
const std::vector<long>& Intensity() const;
```

Description

Get intensity array

1.18.7.3.8 DataPoint::Occluded

```
bool Occluded() const;
```

Description

Get occlusion flag

1.18.7.3.9 DataPoint::Residual

```
const std::vector<long>& Residual() const;
```

Description

Get residual array

1.18.7.3.10 DataPoint::SetOccluded

```
void SetOccluded(bool o);
```

Parameters

| Parameters | Description |
|------------|-------------|
| bool o | flag value |

Description

Set occlusion flag

1.18.7.3.11 DataPoint::Value

```
const std::vector<float>& Value() const;
```

Description

Get value array

1.18.7.4 DataReader

Class Hierarchy

```
TextReader
  DataReader

class DataReader : protected TextReader;
```

File

DataReader.h

Description

Read data from structured text as output by DataWriter (see page 36)

Members**Methods**

| Method | Description |
|---|--|
| ◆ DataReader (see page 35) | Construct for specified stream |
| ◆ LoadAcquisition (see page 35) | Load structured acquisition object as written using DataWriter (see page 36) |
| ◆ LocateObject (see page 35) | Find start of object with specified type and title |
| ◆ ParseFloat32 (see page 36) | Helper to parse string to 32-bit floating point |
| ◆ ParseInt32 (see page 36) | Helper to parse string to 32-bit integer |

Legend

| | |
|---|--------|
| ◆ | Method |
|---|--------|

1.18.7.4.1 DataReader::DataReader

```
DataReader(std::istream& _stream);
```

Parameters

| Parameters | Description |
|-----------------------|---------------------|
| std::istream& _stream | Stream to read from |

Description

Construct for specified stream

1.18.7.4.2 DataReader::LoadAcquisition

```
void LoadAcquisition(DataAcquisition& acquisition, const std::string& title)
throw(TracedException);
```

Parameters

| Parameters | Description |
|------------------------------|---|
| DataAcquisition& acquisition | To fill with acquisition read |
| const std::string& title | The title of the acquisition object to read |

Description

Load structured acquisition object as written using [DataWriter](#) (see page 36)

Exceptions

[TracedException](#) (see page 18) if problems parsing acquisition

1.18.7.4.3 DataReader::LocateObject

```
bool LocateObject(const std::string& objtype, const std::string& title)
throw(TracedException);
```

Parameters

| Parameters | Description |
|----------------------------|-------------------|
| const std::string& objtype | type to look for |
| const std::string& title | title to look for |

Returns

true if found, false if not

Description

Find start of object with specified type and title

Exceptions

[TracedException](#) (see page 18) if there was a problem parsing type or title

1.18.7.4.4 DataReader::ParseFloat32

```
static float ParseFloat32(const std::string& value) throw(TracedException);
```

Parameters

| Parameters | Description |
|--------------------------|-----------------|
| const std::string& value | String to parse |

Returns

parsed value

Description

Helper to parse string to 32-bit floating point

Exceptions

TracedException (see page 18) if could not parse

1.18.7.4.5 DataReader::ParseInt32

```
static long ParseInt32(const std::string& value) throw(TracedException);
```

Parameters

| Parameters | Description |
|--------------------------|-----------------|
| const std::string& value | String to parse |

Returns

parsed value

Description

Helper to parse string to 32-bit integer

Exceptions

TracedException (see page 18) if could not parse

1.18.7.5 DataWriter

Write data to a structured text file

Class Hierarchy

DataWriter

```
class DataWriter;
```

File

DataWriter.h

Description

Writes data in structured text format. This is intended to be a subset of the JSON standard JSON format (<http://www.json.org/>) though some aspects such as escaped character sequences are not supported.

Members

Methods

| Method | Description |
|--------------------------------------|---|
| ◆ AddAcquisitionPacket (see page 37) | Write acquisition packet in structured format - supports 3DResultExt and ADC16 packet layouts |
| ◆ AddVersion (see page 37) | Write version information |
| ◆ BeginAcquisition (see page 37) | Begin an acquisition object. |

| | |
|--------------------------------|--|
| ◆ BeginFile (see page 37) | Write file opening section, just an open bracket |
| ◆ DataWriter (see page 37) | Construct to write to specified stream |
| ◆ EndAcquisition (see page 38) | Write end-of-acquisition characters |
| ◆ EndFile (see page 38) | Complete file to form JSON we can parse |

Legend

| | |
|---|--------|
| ◆ | Method |
|---|--------|

1.18.7.5.1 DataWriter::AddAcquisitionPacket

```
void AddAcquisitionPacket(const codaRTNet::RTNetworkPacket& packet, bool sync_stop = false);
```

Parameters

| Parameters | Description |
|--|--|
| const codaRTNet::RTNetworkPacket& packet | Packet in know format @parma sync_stop Useful for external sync diagnostics - sets an external sync flag in packet |

Description

Write acquisition packet in structured format - supports 3DResultExt and ADC16 packet layouts

1.18.7.5.2 DataWriter::AddVersion

```
void AddVersion(const std::string& title, const codaRTNet::Version& v);
```

Parameters

| Parameters | Description |
|-----------------------------|--|
| const std::string& title | The title to associate with this version |
| const codaRTNet::Version& v | Version using RTNet SDK version class |

Description

Write version information

1.18.7.5.3 DataWriter::BeginAcquisition

Begin an acquisition object.

```
void BeginAcquisition(const std::string& title, const std::map<WORD, float>& device_sample_rate);
```

Parameters

| Parameters | Description |
|--------------------------|--------------------|
| const std::string& title | The title to write |
| Sample | rate list to write |

Description

Will write a type and title, followed by sample rates for each device, followed by start of packet list.

1.18.7.5.4 DataWriter::BeginFile

```
void BeginFile();
```

Description

Write file opening section, just an open bracket

1.18.7.5.5 DataWriter::DataWriter

```
DataWriter(std::ostream& _stream);
```

Parameters

| Parameters | Description |
|--|--------------------|
| <code>std::ostream& _stream</code> | Stream to write to |

Description

Construct to write to specified stream

1.18.7.5.6 DataWriter::EndAcquisition

```
void EndAcquisition();
```

Description

Write end-of-acquisition characters

1.18.7.5.7 DataWriter::EndFile

```
void EndFile();
```

Description

Complete file to form JSON we can parse

1.18.7.6 TextReader

Read structured text files

Class Hierarchy

TextReader

```
class TextReader;
```







File

TextReader.h


Description

This is designed to read parts of the structured text files which are output from DataWriter (see page 36). They use a subset of the JSON format (<http://www.json.org/>).

Members**Methods**

| Method | Description |
|---|---|
|  NextNonSpace (see page 38) | Skip to next non-space and read it |
|  ParseKeyValue (see page 39) | Parse a key and value separated by colon and arbitrary numbers of space characters |
|  ParseQuotedString (see page 39) | Read from quote character until and including the next quote character and extract the string in between. |
|  ParseValue (see page 39) | Parse a value. |
|  SkipSpace (see page 40) | Skip space |
|  TextReader (see page 40) | Construct to read from given stream |

Legend

| | |
|---|--------|
|  | Method |
|---|--------|

1.18.7.6.1 TextReader::NextNonSpace

```
int NextNonSpace();
```

Returns

The value of next non-space character

Description

Skip to next non-space and read it

1.18.7.6.2 TextReader::ParseKeyValue

Parse a key and value separated by colon and arbitrary numbers of space characters

```
void ParseKeyValue(std::string& key, std::string& value) throw(TracedException);
```

Parameters

| Parameters | Description |
|--------------------|----------------------|
| std::string& key | Filled with key read |
| std::string& value | Filled with value |

Description

Parse a key-value pair. Supports quoted and non-quoted values. For example:

"mykey" : "myvalue"

Or:

"some number" : 3

Exceptions

TracedException (see page 18) if stream is not positioned before such a pair

1.18.7.6.3 TextReader::ParseQuotedString

```
void ParseQuotedString(std::string& str) throw(TracedException);
```

Parameters

| Parameters | Description |
|------------------|---|
| std::string& str | Filled with the string between the quotes |

Description

Read from quote character until and including the next quote character and extract the string in between.

Exceptions

TracedException (see page 18) if the stream is not positioned just before a quote

1.18.7.6.4 TextReader::ParseValue

Parse a value.

```
void ParseValue(std::string& value) throw(TracedException);
```

Parameters

| Parameters | Description |
|------------|----------------------------|
| str | Filled with the value read |

Description

If the next non-space character is a quote, this will read up to and including the next quote and return the string in between (like ParseQuotedString (see page 39)). But if the next non-space character is not a quote then it will read up until any space or control character is found, and so can read numbers which are not generally enclosed in quotes

Exceptions

TracedException (see page 18)

1.18.7.6.5 TextReader::SkipSpace

```
void SkipSpace();
```

Description

Skip space

1.18.7.6.6 TextReader::TextReader

```
TextReader(std::istream& _stream);
```

Parameters

| Parameters | Description |
|-----------------------|----------------|
| std::istream& _stream | Stream to read |

Description

Construct to read from given stream

1.18.7.7 CommandListReader

Class Hierarchy

```
TextReader
  CommandListReader

class CommandListReader : protected TextReader;
```

File

CommandListReader.h

Description

Read (see page 41) list of commands from structured text file

Members

Methods

| Method | Description |
|---------------------------------|---|
| CommandListReader (see page 40) | Construct to read from specified stream |
| Read (see page 41) | Parse according to registered command factories and build command instances |

Legend

| | |
|---|--------|
|  | Method |
|---|--------|

1.18.7.7.1 CommandListReader::CommandListReader

```
CommandListReader(std::istream& _stream);
```

Parameters

| Parameters | Description |
|-----------------------|----------------|
| std::istream& _stream | Stream to read |

Description

Construct to read from specified stream

1.18.7.7.2 CommandListReader::Read

```
void Read(CommandList& commandlist) throw(TracedException);
```

Parameters

| Parameters | Description |
|--------------------------|--|
| CommandList& commandlist | Command (see page 20) list to populate |

Description

Parse according to registered command factories and build command instances

Exceptions

TracedException (see page 18) if structured text format was not as expected

1.18.8 Tests

These are automated self-tests for certain sections of code.

Functions

| Function | Description |
|-----------------------------------|---|
| SelfTest (see page 41) | Run all self-tests |
| SelfTestCommandLine (see page 42) | Test command line parsing |
| SelfTestDataReader (see page 42) | Test data reader |
| SelfTestDataWriter (see page 42) | Test data writer |
| SelfTestScripts (see page 42) | Test command script parsing and execution |

Macros

| Macro | Description |
|--------------------------------|---|
| SELF_TEST_ASSERT (see page 43) | Self-test failure if condition is false. The code given for the condition will be logged as the failure message |
| SELF_TEST_FAIL (see page 43) | Indicate self-test failure with specified message |

Variables

| Variable | Description |
|----------------------------------|--|
| SelfTestModuleName (see page 43) | Identifier to use for module name if an exception is thrown during self-test |

1.18.8.1 SelfTest

```
void SelfTest(ResultLog& results) throw(TracedException);
```

File

SelfTest.h

Parameters

| Parameters | Description |
|--------------------|-------------------------------|
| ResultLog& results | Log list of tests to this log |

Description

Run all self-tests

Exceptions

TracedException (see page 18)

1.18.8.2 SelfTestCommandLine

```
void SelfTestCommandLine() throw(TracedException);
```

File

SelfTestCommandLine.h

Description

Test command line parsing

Exceptions

TracedException (see page 18)

1.18.8.3 SelfTestDataReader

```
void SelfTestDataReader() throw(TracedException);
```

File

SelfTestDataReader.h

Description

Test data reader

Exceptions

TracedException (see page 18)

1.18.8.4 SelfTestDataWriter

```
void SelfTestDataWriter() throw(TracedException);
```

File

SelfTestDataWriter.h

Description

Test data writer

Exceptions

TracedException (see page 18)

1.18.8.5 SelfTestScripts

```
void SelfTestScripts() throw(TracedException);
```

File

SelfTestScripts.h

Description

Test command script parsing and execution

Exceptions

TracedException (see page 18)

1.18.8.6 SelfTestModuleName

```
const char* SelfTestModuleName;
```

File

SelfTestMacros.h

Description

Identifier to use for module name if an exception is thrown during self-test

1.18.8.7 SELF_TEST_ASSERT

```
#define SELF_TEST_ASSERT(condition) { if (!(condition)) { SELF_TEST_FAIL(#condition); } }
```

File

SelfTestMacros.h

Parameters

| Parameters | Description |
|------------|-------------------|
| condition | Condition to test |

Description

Self-test failure if condition is false. The code given for the condition will be logged as the failure message

1.18.8.8 SELF_TEST_FAIL

```
#define SELF_TEST_FAIL(_message) STOP(SelfTestModuleName, _message)
```

File

SelfTestMacros.h

Parameters

| Parameters | Description |
|------------|--------------------|
| _message | The message to log |

Description

Indicate self-test failure with specified message

Index

A

- AcqContinuous 2
 - AcqContinuous 3
 - RunRTNet 3
- AcqContinuousBuffered 4
 - AcqContinuousBuffered 4
 - RunRTNet 4
- AcqContinuousBufferedFactory 5
 - AcqContinuousBufferedFactory 5
- AcqContinuousFactory 3
 - AcqContinuousFactory 3
- AcqSingleShot 5
 - AcqSingleShot 6
 - RunRTNet 6
- AcqSingleShotFactory 6
 - AcqSingleShotFactory 7
- Alignment 1
 - Alignment 1
 - RunRTNet 1
- AlignmentFactory 2
 - AlignmentFactory 2

C

- CodaMode 7
 - CodaMode 7
 - RunRTNet 7
- CodaModeFactory 8
 - CodaModeFactory 8
- Command 20
 - ~Command 21
 - class CommandFactory 23
 - Command 21
 - Comment 21
 - Comment1 21, 22
 - GetParameter 22
 - HaveRequiredParameters 22
 - Name 22
 - RegisterParameter 22
 - Run 22

- COMMAND_STOP 28
- COMMAND_STOP2 28
- CommandFactory 23
 - Build 23
 - CommandFactory 23
- CommandFactoryRTNet 12
 - Build 13
 - CommandFactoryRTNet 13
- CommandFactoryRTNetData 14
 - Build 15
 - CommandFactoryRTNetData 15
- CommandLineOptions 15
 - CommandLineOptions 16
 - OptionRegister 16
 - Parse 16
 - RegisterOption 16
- CommandList 23
 - ~CommandList 24
 - BuildCommand 24
 - CommandList 24
 - RegisterCommand 24
 - RunAll 25
- CommandListReader 40
 - CommandListReader 40
 - Read 41
- CommandRTNet 11
 - class CommandFactoryRTNet 12
 - Client 12
 - Run 12
 - RunRTNet 12
- CommandRTNetData 13
 - class CommandFactoryRTNetData 14
 - DataOutput 14
- Commands 20

D

- DataAcquisition 29
 - AddPacket 29
 - DataAcquisition 30
 - GetDeviceSampleRate 30
 - Packets 30
 - SetDeviceSampleRate 30

| | |
|-------------------------|----------------------|
| DataPacket 30 | ~Parameter 25 |
| = 31 | IsSet 25 |
| AddPoint 31 | Parameter 25 |
| DataPacket 32 | Parse 26 |
| device 31 | ParameterInteger 26 |
| page 31 | ~ParameterInteger 26 |
| Points 32 | ParameterInteger 26 |
| tick 31 | Parse 27 |
| DataPoint 32 | Value 27 |
| = 33 | ParameterString 27 |
| AddIntensity 33 | ~ParameterString 27 |
| AddResidual 33 | ParameterString 27 |
| AddValue 33 | Parse 28 |
| DataPoint 33 | Value 28 |
| Intensity 34 | |
| Occluded 34 | |
| Residual 34 | |
| SetOccluded 34 | |
| Value 34 | |
| DataReader 34 | |
| DataReader 35 | |
| LoadAcquisition 35 | |
| LocateObject 35 | |
| ParseFloat32 36 | |
| ParseInt32 36 | |
| DataWriter 36 | |
| AddAcquisitionPacket 37 | |
| AddVersion 37 | |
| BeginAcquisition 37 | |
| BeginFile 37 | |
| DataWriter 37 | |
| EndAcquisition 38 | |
| EndFile 38 | |

F

Framework 15

I

IO 29

P

Parameter 25

R

ResultLog 17

- Log1 17
- Log2 17
- LogException 18
- ResultLog 18

RTNetDemoOptions 10

- CommandFile 10
- DataFile 10
- Parse 10
- RTNetDemoOptions 11
- Server 11

S

SELF_TEST_ASSERT 43

SELF_TEST_FAIL 43

SelfTest 41

SelfTestCommandLine 42

SelfTestDataReader 42

SelfTestDataWriter 42

SelfTestModuleName 43

SelfTestScripts 42

StartSystem 8

- RunRTNet 9
- StartSystem 9

StartSystemFactory 9

- StartSystemFactory 9

STOP 19

STOP2 19

T

Tests 41

 StreamReader 38

 NextNonSpace 38

 ParseKeyValue 39

 ParseQuotedString 39

 ParseValue 39

 SkipSpace 40

 StreamReader 40

 TracedException 18

 Report 19

 TracedException 19